

.

From: support@agilent.com
Sent: Wednesday, October 10, 2012 11:44 AM
To: Hanchett, James (DPH)
Subject: MS ChemStation Error message from George, SR# 1-4693852882

Hi Jim,

Regarding the error "File 51 not open for OUTPUT or APPEND " you were seeing earlier this week, did deleting the file QA.LOG fix the problem?

Regards,

George Picarello
Customer Service Engineer
Online Technical Support Team
Agilent Technologies, Inc.
Phone: 1-800-227-9770
www.home.agilent.com



This email may contain confidential information and is meant for the intended recipient(s) only. If you are not an intended recipient or someone authorized to receive e-mail on behalf of an intended recipient, please reply to the sender that the email was misdirected and delete the email.